



GALIX

BIOMEDICAL INSTRUMENTATION, INC.

GBI-3SM: Digital ECG Holter Recorder

- 24/48/72 hours of ECG recording with only one AAA alkaline battery
- Very High Resolution ECG
- Small, Ultra Light, Ultra Slim design
- Pacemaker Detection
- Graphic ECG Display



Other Features:

- Dimensions: 3"W x 2.2"H x 0.65"D (76 x 57 x 16 mm)
- Weight: approximately 2 ounces with battery (60 grams)
- Graphic display for all 3 ECG channels and data
- Programmable Sampling Rate: 256 - 1000 samples/sec
- Pacemaker pulse detection: as small as 0.25 ms width
- 24, 48, 72 hours with one AAA alkaline battery
- Frequency response: 0.05 – 100 Hz.
- Recording Media: Secure Digital Memory Card
- Event Button with tactile and audible feedback
- Versatile, user friendly interactive screen with ECG display
- Environmental directives: RoHS compliant