



GALIX

BIOMEDICAL INSTRUMENTATION, INC.

PHYSIO-28

Electrophysiology Polygraph with built-in Stimulator

- Low cost, yet powerful Electrophysiology System
- Very High Resolution ECG Data Acquisition
- Digital acquisition of 12 surface ECG and 16 intracardiac channels
- Extremely reliable and versatile built-in Stimulator, with stimulation channel selected by software.



Other Features:

- Up to three independent extrastimuli with automatic coupling delay increment / decrement and multiprogrammable pacing
- Continuous data storage on recordable media
- Studies can be stored completely or partially. These studies can be retrieved at any time for comparison and/or further analysis purposes
- Can be installed on a standard PC or a Notebook

2555 Collins Ave., Suite C-5, Miami Beach, FL 33140, USA. Phone: 305-534-5905. Fax: 305-534-8222

sales@galix-gbi.com

www.galix-gbi.com